

Mounting of the mattress


RotoBed®Home comes with a mattress consisting of a main part and a foot part.

The main mattress has 5 buttons on the back. First place the three lower buttons in the grooves at the bottom of the bed. Slide the mattress towards the foot end. Now press the two buttons in the centre of the mattress into the corresponding holes in the middle of the bed.

Similarly, the foot mattress has two buttons, which fit into the two round holes in the bed.

Mounting of side rails

RotoBed®Home can be delivered with side rails made of metal or oak laminate. Regardless of the material the side rails are mounted and removed in the same way:



Press the small fitting underneath while lifting the side rail.

The side rails can simply be put back into place and will automatically lock.

Do not leave the side rails upright and loose in the bracket. This may cause squeeze risk and damage to the bed during rotation.

Note: RotoBed®Home should not be used by people without an introduction in the use of the bed.

RotoBed®Home is CE-marked according to directive (EU) 2006/42/EF.

RotoBed®Home should be operated exclusively with the remote control. Be sure that children are aware of the squeeze risk when adjusting the bed's positions. When the bed's leg part is raised, do not overload it with heavy objects or persons.

Quick guide for RotoBed® Home



Entry: Get seated as close to the backrest as possible for the best comfort. Hold down the **Sleeping Position** key until the bed is rotated and horizontal. Now you can adjust backrest and leg part as wished.




Exit: Hold down the **Chair Position** key until the bed stops in an exit position. To facilitate exit you can, by means of the remote control, move the backrest the last bit forward.



Note: As standard, RotoBed®Home is programmed to run to an oblique exit position to facilitate transfer to a wheelchair or a walker. If you wish a perpendicular exit position, please see user manual.

Remote control



RotoBed®Home can turn either right or left. Therefore, the remote control can also be connected to both the right and the left side of the bed. When the remote control is not in use, place it on the side rail.